



1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23

Parenteral and Ophthalmic Drug Products (PODP) Leachables and Extractables Working Group

**Issued and Effective
September, 2011**

Study Protocol – Stage 2

Experimental Protocol for Simulation Study of Blow-Fill-Seal (BFS) PODP Container Closure Systems

<u>TABLE OF CONTENTS</u>	
24	
25	
26	I. INTRODUCTION 3
27	
28	II. PURPOSE AND SCOPE OF WORK (STUDY PROTOCOL STAGE 2) 5
29	
30	III. REGULATORY STATUS 6
31	
32	IV. SAFETY AND ENVIRONMENTAL IMPACT 6
33	
34	V. TEST SYSTEM 7
35	
36	VI. CHEMICALS AND EQUIPMENT 7
37	
38	A. Extraction Solvents 8
39	B. Additional Chemicals 8
40	C. Extraction Equipment 8
41	D. Analytical Instrumentation 8
42	
43	VII. EXTRACTION PROCEDURES 9
44	
45	A. Simulating Extraction Solvents 9
46	B. Accelerated Extraction Conditions 9
47	C. Exaggerating Factors 9
48	
49	VIII. ANALYTICAL METHODS 10
50	
51	A. General 10
52	B. System Suitability 10
53	B. Gas Chromatography (GC) 18
54	1. General 18
55	2. Sample Preparation 18
56	3. Operating Conditions 22
57	4. General Comments 22
58	C. High Performance Liquid Chromatography (HPLC) 22
59	1. General 22
60	2. Sample Preparation 22
61	3. Operating Conditions 22
62	4. General Comments 23
63	E. Headspace GC/MS Analysis 23
64	1. General 23
65	2. Sample Preparation 23
66	3. Operating Conditions 24
67	4. General Comments 24

68	F.	Inductively Coupled Plasma Atomic Spectroscopy (ICPAS)	25
69	1.	General	25
70	2.	Sample Preparation	25
71	3.	Operating Conditions	25
72	4.	General Comments	26
73			
74	IX.	DATA EVALUATION AND REPORTING	26
75			
76	A.	Qualitative Analysis	26
77	B.	Semi-Quantitative Measurement	26
78			
79	X.	GLOSSARY	27
80			
81	XI.	REFERENCES	27
82			
83			

84 **I. Introduction**
85

86 It has been well established that substances extracted by drug products from their container
87 closure systems can affect the drug product's safety and efficacy. Regulatory guidance has
88 provided some recommendations regarding the analysis and toxicological safety assessment (*i.e.*,
89 qualification) of such substances. Thus, for example, the FDA issued *Container Closure Systems*
90 *for Packaging Human Drugs and Biologics – Chemistry, Manufacturing and Controls (CMC)*
91 *documentation Guidance for Industry* in May 1999¹. In addition, the European Medicines
92 Agency (EMA) issued its *Guideline on Plastic Immediate Packaging Materials* in May 2005.²
93 Specific Guidance for Orally Inhaled and Nasal Drug Products (OINDP) is contained in two
94 CMC Guidances addressing OINDP¹: (i) the draft *Guidance for Industry, Metered Dose Inhaler*
95 *(MDI) and Dry Powder Inhaler (DPI) Drug Products, Chemistry, Manufacturing, and Controls*
96 *Documentation* (November, 1998); and (ii) the *Guidance for Industry, Nasal Spray and*
97 *Inhalation Solution, Suspension, and Spray Drug Products Chemistry, Manufacturing, and*
98 *Controls Documentation* (July, 2002).
99


100 In September 2006, the Product Quality Research Institute (PQRI) issued a Recommendation
101 entitled "Safety Thresholds and Best Practices for Extractables and Leachables in Orally Inhaled
102 and Nasal Drug Products"³. This Recommendation provided a scientific rationale and process to
103 identify, quantify and establish the biological safety (*i.e.* qualify) of leachables and/or
104 extractables where appropriate, in OINDP. Included in this Recommendation were experimental
105 protocols, and the results thereof, for establishing Best Demonstrated Practices for the
106 performance of Controlled Extraction Studies, specifically relevant of the OINDP dosage forms.
107

108 The PQRI Parenteral and Ophthalmic Drug Products (PODP) Leachables and Extractables
109 Working Group is developing, executing and reporting experimental studies as the means of
110 establishing Best Demonstrated Practices for the performance of Chemical Assessments
111 specifically relevant for PODP container closure systems and dosage forms. Figure 1, The
112 Chemical Assessment Triad, illustrates the Chemical Assessment Process. The PODP Stage 1
113 study considered the process of Material Characterization; specifically the processes by which a
114 Controlled Extract is generated, by which a Controlled Extract is analyzed and by which the test
115 results are evaluated and interpreted. This Stage 2 study considers the process of performing a
116 Simulation Study, specifically establishing the extractables profile of an experimental container
117 closure system constructed from some of the materials that were characterized in the Stage 1
118 study. This experimental container closure system specifically mimics a Blow-Fill-Seal (BFS)
119 packaging system, such as those used with many ophthalmic products, consisting of a BFS
120 bottle, its associated cap, a closure gasket and an affixed printed label.
121

122 This experimental protocol will be used by all participating laboratories and investigators.
123
124
125

¹ Available at <http://www.fda.gov/downloads/Drugs/GuidanceComplianceRegulatoryInformation/Guidances/ucm070551.pdf>
² Available at <http://www.emea.europa.eu/pdfs/human/qwp/435903en.pdf>
³ Available at http://pqri.org/pdfs/LE_Recommendations_to_FDA_09-29-06.pdf

126 **Figure 1.** Generalized Chemical Assessment, The Chemical Assessment Triad. A simulation
127 study may be an appropriate and effective bridge between the material
128 characterization process, which establishes extractables that are tentative leachables
129 and product assessment, which measured confirmed leachables. As the name
130 suggests, a simulation study seeks to mimic the product assessment by using
131 simulating solvents to facilitate the analytical tasks and using extraction conditions
132 which accelerate the product contact conditions. The simulation study may be the
133 basis of a preliminary toxicological assessment and can be used to establish target
134 leachables to measure during product assessment.



135
136
137
138
139
140
141
142
143
144
145
146
147
148
149
150
151
152
153
154
155
156
157
158
159
160
161
162
163
164
165
166
167
168
169
170
171
172
173
174
175
176
177

178 II. Purpose and Scope of Work (Study Protocol Stage 2)

179
180 The purpose of the experiments outlined in this protocol is to generate data from a Simulation
181 Studies, which the Working Group will use to investigate its hypotheses:
182

- 183 1. Threshold concepts that have been developed for safety qualification of leachables in
184 OINDP can be extrapolated to the evaluation and safety qualification of leachables in
185 PODP, with consideration of factors and parameters such as dose, duration, patient
186 population and product dependent characteristics unique to various PODP types.
187
- 188 2. The science-based best demonstrated practices established for the OINDP pharmaceutical
189 development process can be extrapolated to PODP container closure systems.
190
- 191 3. Threshold and best practices concepts can be integrated into a comprehensive process for
192 characterizing container closure systems with respect to leachable substances and their
193 associated impact on PODP safety.
194

195 The Simulation Study will be performed following the general methodologies contained in this
196 protocol. The Test System, designed to mimic a Blow-Fill-Seal (BFS) packaging system, will be
197 filled with various simulating extraction solvents which are intended to mimic aqueous PODP
198 drug products. The Test System will be exposed to accelerated storage conditions that mimic a
199 post-filling shelf life of 2 years at ambient temperature. Consistent with the BFS process, the
200 filled Test System will not be exposed to terminal or auxiliary sterilization. At certain times
201 during accelerated storage, the simulating extraction solvents will be characterized to reveal
202 organic and inorganic extractables, thus establishing the Test System's Extractables Profile. The
203 intent of this Stage 2 assessment is to generate the fundamental information from which Best
204 Demonstrated Practices can be derived; it is not the intent of this Stage 2 assessment to
205 prospectively establish the practices used in this study as the Best Demonstrated Practices
206 themselves. It can be reasonably expected that this Extractables Profile will include tentative
207 extractables revealed in the Stage 1, Material Characterization, study,
208

209 As no single analytical technique can be used to identify and quantify all extractables, a variety
210 of methods will be utilized in this protocol to maximize the likelihood that all predominant
211 extractable compounds associated with the test articles are accounted for and appropriately
212 evaluated. Overlap between methods will supply corroborating data that demonstrate the validity
213 of the procedures. To provide a full analytical survey of possible analytes the following strategy
214 will be employed:
215

- 216 1. Gas Chromatography with appropriate sampling/injection and detection strategies
217 e.g. Flame Ionization Detection (GC/FID) and Mass Spectrometry (GC/MS)] for
218 identification and assessment of volatile and semi-volatile extractables.
- 219 2. High Performance Liquid Chromatography with appropriate detection strategies
220 [e.g. Diode Array Detection (HPLC/DAD), Mass Spectrometry (LC/MS)] for
221 identification and assessment of relatively polar and non-volatile extractables.

- 222 3. Gas chromatography with headspace sampling (HS-GC) for identification and
223 assessment of volatile extractables.
224 4. Inductively Coupled Plasma/Mass Spectrometry (ICP/MS) and/or Inductively
225 Coupled Plasma/Atomic Emission Spectroscopy (ICP/AES) to detect single
226 elements in the extracts (i.e. metals).
227

228 While analytical tests and measurements, such as pH, UV absorbance, and total organic carbon
229 (TOC), can provide insight into the general chemical nature and amount of extracted substances,
230 they do not directly provide information for the identification and/or quantitation of individual
231 extractables and thus will not be utilized in this study.
232

233 Studies designed to assess recovery (i.e. mass balance) for individual extractables relative to the
234 known formulations of chemical additives in the various test articles, or reproducibility of
235 extractables profiles for multiple “batches” of any particular test article are not within the scope
236 of this Stage of the test protocol. Additionally, the extraction procedures, analytical
237 techniques/methods, and analysis conditions described in this experimental test protocol will not
238 be fully and rigorously validated. Nevertheless, the scientific credibility of the data generated in
239 this study shall be established via the utilization of system suitability testing with all the analysis
240 methods and by the expert review of the generated data. Finally, “special case” classes of
241 extractables that have defined and highly specific analytical methods that are generally accepted
242 and commonly used for their identification and quantitative assessment will not be considered in
243 this study.
244

245 **III. REGULATORY STATUS**

246
247 This experimental test protocol will be conducted in the spirit of Good Laboratory Practices and
248 Good Manufacturing Practices (GXP) requirements. All experiments shall be documented based
249 on the appropriate GXP compliance systems in a participating laboratory. Any changes or
250 clarifications that a participating laboratory makes to this test protocol shall be documented as
251 appropriate, and discussed/approved by the Study Coordination as appropriate.
252

253 **IV. SAFETY AND ENVIRONMENTAL IMPACT**

254
255 Chemicals and reagents used in this study (e.g. organic solvents commonly used to enhance
256 solubility of lipophilic targets and to increase transport of small molecules out of complex
257 matrices) may be flammable and/or pose short-term and long-term environmental health risks.
258 Care must be exercised with their use. Consult the Material Safety and Data Sheet (MSDS) for
259 appropriate personal protection and disposal. Safety risks associated with the various processes
260 and procedures performed in this study may exist and should be understood and managed using
261 such strategies as environmental control and personal protection.
262

263 **V. TEST SYSTEM**

264
265 The Test System is designed to mimic a Blow-Fill-Seal (BFS) packaging system, such as that used
266 for ophthalmics and small volume parenterals (SVP). In general such a system consists of the

267 BFS container (in this case a bottle), one or more closures (in this case the bottle's associated cap
 268 and a gasket/liner) and some type of labeling (in this case a printed label). Due to procurement
 269 limitations and confidentiality issues, a Test System that truly represents commercial BFS
 270 packaging systems was not obtained. Rather, a test system was loosely constructed from
 271 materials which had been characterized in the Stage 1 study. Specifically, the Test System
 272 consists of:

- 273
- 274 1. A low density polyethylene (LDPE) bottle, with a PP screw cap,
- 275 2. A rubber gasket which is used as a liner inside the plastic cap, and
- 276 3. A printed adhesive label, which is affixed to the bottle's outer surface.
- 277

278 Identifying information related to these Test System Components is contained in Table 1.
 279

TABLE 1. TEST ARTICLE.			
MATERIAL TYPE	MATERIAL APPLICATION	MATERIAL FORMAT	DESCRIPTION
Low density polyethylene (LDPE)	Bottle/ Vial	Bottle	4 oz LDPE, part B347A (Container & Packaging Supply)
Polypropylene (PP)	Cap	Cap	PP, Part L764(Container & Packaging Supply)
Adhesive Label	Label on Container Surface	Label Sheets	<u>Substrate</u> : Unknown <u>Adhesive</u> : Acrylic polymer(s), residual monomers, water, ammonia (99.55%); wetting agent, Surfynol 336, at 0.4% containing CAS 577-11-7 (> 25%), CAS 9014-85-1 (> 25%); Biocide, Kathon LX, at 0.05% containing Chloro-2-methyl-4-isothiazolin-3-one (CAS 26172-55-4), 1.1-1.4%,2-Methyl-4-isothiazolin-3-one (CAS 2682-20-4), 0.3 - 0.5%, Magnesium Chloride (CAS 7786-30-3), 1.0 - 1.2%, Magnesium nitrate (CAS 10377-60-3), 1.4 - 2.0% Copper nitrate (CAS 3251-23-8) 1,500 - 1,700 ppm, Water, 95 - 97% <u>Printing ink</u> : Irgacure 369 (CAS 119313-12-1) and Irgacure 1173 (CAS 7473-98-5), photoinitiators; Trimethylolpropane triacrylate (TMPTA, CAS 15625-89-5), Tripropylene glycol diacrylate (TPGDA, CAS 42978-66-5), Glycerol propoxy triacrylate (GPTA, CAS 52408-84-1), monomers; HQME/Mequinol (CAS 150-76-5), stabilizer; Carbon black (CAS 1333-86-4),Phthalo blue (CAS 147-14-8),Carbazole violet (CAS 215247-95-3), pigments <u>Varnish</u> : Unknown
Rubber (Elastomer) (RE)	Closures	Gasket/liner	Brominated isobutylene isoprene copolymer (57.3%); calcined aluminum silicate, 38.2%, titanium dioxide, 1.2%; paraffinic oil, 1.2%; zinc oxide, 0.6%; polyethylene, 0.6%; SRF Carbon block mixture, 0.4%; calcined magnesium oxide, 0.3%; 4,4'-dithiodi-morpholine/polyisobutylene, 0.3%

280
 281 **VI. CHEMICALS AND EQUIPMENT**

282
 283 Extraction and analytical methods were chosen and designed to utilize chemicals, apparatus, and
 284 instrumentation available in typical laboratories routinely involved with this type of study.

285 **A. Extraction Solvents**

286
287 Chemicals required for the use as, or preparation of, extraction solvents, are as follows:
288

- 289 • Laboratory research grade water or Water for Injection (WFI), appropriately sourced,
290 collected and stored to minimize background levels of extraneous substances.
- 291 • Potassium chloride
- 292 • Hydrochloric acid, 0.1 N
- 293 • Sodium phosphate monobasic
- 294 • Sodium phosphate dibasic
- 295 • Sodium hydroxide, 1 N
- 296 • Isopropyl alcohol (glass bottled; IPA)
- 297 • pH calibration buffers; pH 1.68, 4.01, 9.18 and 12.48 (saturated calcium hydroxide)

298
299 The preparation of several of these extraction solvents is as follows:
300

- 301 • Water at pH 2.5 (HCl/KCl mixture): The KCl solution is prepared at 0.01M. Weigh
302 1.5 grams of KCl into a 2.0 L vol flask containing 1500 mL water. Add 60 mL 0.1 N
303 HCl. Dilute to volume with water. This final solution is 0.01 M KCl and 0.003 M
304 HCl, which should have a pH of 2.5.
- 305
- 306 • Water at pH 9.5: Weigh 1.24 grams sodium phosphate monobasic and 18.7 grams of
307 sodium phosphate dibasic, transfer to an appropriate vessel, and dissolve in 2 liters of
308 water. The pH of this solution is reported to be 8.0 Titrate with 1 N NaOH to get a
309 pH of 9.5. This solution is 0.0045 M monobasic and 0.066 M dibasic.
- 310
- 311 • IPA/Water (1/1): Mix equal volumes of IPA and water.

312
313 **B. Additional Chemicals**

- 314
- 315 • Analytical reagents required to perform the analytical testing.
- 316 • Reference and/or Internal standards required to perform the analytical testing.

317
318 **C. Extraction Equipment**

- 319
- 320 • Oven with operating range of 30 to 50 °C; explosion proof

321
322 **D. Analytical Instrumentation**

- 323
- 324 • Gas chromatograph equipped with a Flame Ionization Detector (GC/FID)
- 325 • Gas chromatograph equipped with a Mass Spectrometer (GC/MS). GC systems that employ
326 flow splitting to accomplish FID and MS detection in tandem could be used in this study.
- 327 • Headspace Sampler/Injector (HS) for GC/MS Instrumentation.
- 328 • Liquid chromatograph equipped with a photodiode array detector

- 329 • Liquid chromatograph equipped with an APCI (Atmospheric Pressure Chemical Ionization)
330 capable Mass Spectrometer (LC/MS). Preference is given to LC systems that are capable of
331 both DAD and MS detection. Additional detectors (e.g. corona assisted discharge detectors,
332 evaporative light scattering) may be used as appropriate.
333 • Inductively Coupled Plasma Mass Spectrometer (ICP/MS, preferred) and/or Inductively
334 Coupled Plasma Atomic Emission Spectrometer (ICP/AES)
335

336 VII. EXTRACTION PROCEDURE

337 A. Simulating Extraction Solvents

338 In the PQRI OINDP studies, extractions were performed on each test article using three solvents
341 representing a range of polarity, specifically

- 342
- 343 • methylene chloride (dichloromethane)
 - 344 • 2-propanol (isopropanol, IPA)
 - 345 • hexane (n-hexane, not hexanes).
- 346

347 This was appropriate in the case of OINDP given the nature of the drug vehicles used in those
348 types of products (organic solvents) and the conditions of contact between the drug vehicles and
349 the container closure system (continuous direct contact over shelf life).

350 While the use of such extraction solvents may be relevant for PODP products, a significant
351 portion of PODP products are water-based and the three solvents previously employed do not
352 address the unique solubilizing properties of water and aqueous buffer systems. Thus in the case
353 of PODP, the OINDP solvents will be replaced by aqueous extraction media. These aqueous
354 extraction media, and their associated justification, include

- 355
- 357 * Water at pH 2.5 (HCl/KCl mixture); justification, few therapeutic products are lower than pH
358 2.5.
 - 359 * Water at pH 9.5 (Phosphate buffer); justification, few therapeutic products are higher in pH
360 than 9.5.
 - 361 * 1/1 IPA/water; justification; simulates aqueous formulations containing solubilizing agents.
 - 362 * Water alone.
- 363

364 B. Accelerated Extraction Conditions

365

366 It is generally well-established that storage at 40°C for 6 months accelerates a product shelf life
367 of 2 years at ambient temperature. Such acceleration will be utilized in this study. In order to
368 allow for trending of the extractables profile versus time, this study will utilize three test
369 intervals; after 1, 2 and 6 months of storage. Thus individual test systems will be filled, stored
370 for either 1, 2 or 6 months and then profiled for extracted substances. The filling will be
371 staggered so that all the filled Test Articles “mature” at the same time, facilitating the analysis of
372 the fill solutions. This will be accomplished by filling the test articles at three different times.
373 One group of test articles will be filled at the initiation of the study and placed into storage.

374 These would be the 6 month samples. After four months of storage has occurred, a second group
375 of test articles will be filled and placed into storage. These would be the 2 months samples.
376 After five months of storage, a third group of samples will be filled and put into storage. These
377 would be the 1 month samples. After 6 months of storage, all units would be removed from
378 storage and tested.

379

380 C. Exaggerating Factors

381

382 The weights and surface areas of the Test System components will be measured so that the Test
383 System can be compared to commercial BFS packaging systems either on the basis of mass or
384 surface area.

385

386 VIII. ANALYTICAL METHODS

387

388 A. General

389

390 The analytical screening methods that will be utilized to discover, identify and quantitate
391 extracted substances will be similar to the same analytical methods that we used to provide the
392 data in the Stage 1 study. As considerable experience was gained with these methods during the
393 Phase 1 study, the actual operating conditions used to support this Phase 2 study may be
394 somewhat different than the conditions enumerated in this Protocol. Such differences will be
395 noted in the Final Report associated with this study.

396 It is never-the-less expected that all analyses performed in this study will be performed with
397 systems and operating parameters that meet the system suitability requirements given herein.

398

399 It is the general objective that all organic extractables that are present in the simulating extracts
400 at a concentration of 0.1 µg/mL (ppm) will be confidently identified and effectively quantitated.
401 It is the general objective that all extracted metals and trace elements that are present in the
402 simulating extracts at a concentration of 0.01 µg/mL (ppm) will be confidently identified and
403 effectively quantitated. It is noted that testing of extraction blanks (portions of simulating
404 solvent stored in inert vessels) allows one to differentiate between analytical artifacts (which are
405 present in both the blanks and the extracts at roughly similar levels) and extractables, which are
406 present in the extracts at levels significantly higher than they are in the blanks.

407

408 B. System Suitability

409

410 All testing performed in support of this Protocol shall include appropriate system suitability
411 assessment. Demonstration of system suitability will be accomplished according to the following
412 three-step approach:

413

414 Step 1: Each participating laboratory will ensure that analytical instrumentation is in proper
415 condition and will demonstrate instrument suitability by following its proprietary (in-
416 house) procedures.

417

418 Step 2: Each participating laboratory will follow the procedures defined in this Protocol which
419 involve the characterization of specified test mixtures by GC, HS-GC, LC and ICP.
420 The test mixtures are suitable to demonstrate adequate and effective analytical
421 performance (for example, separation efficiency, selectivity and sensitivity). All
422 generated system suitability data will be evaluated with regard to the required
423 specifications/acceptance criteria.

424
425 Step 3: Internal Standardization. Specifically for the GC methodology, the extracts will be
426 supplemented by introducing a surrogate internal standard and an injection standard.
427 Analysis of these standards complements system suitability testing by providing a
428 means of establishing the effectiveness of sample preparation/sample introduction
429 processes. The use of internal standards is discussed in the section describing the actual
430 GC analysis of the extracts.

431
432 Table 2 presents a list of system suitability analytes for GC and HPLC based analytical
433 techniques.

434
435 System suitability testing for the ICP trace element analysis shall include the preparation and
436 testing of a system suitability test mixtures that contains all the targeted elements listed
437 previously at a concentration of 0.25 µg/ml. System suitability testing shall consist of the
438 demonstration that all elements can be detected at the prepared concentration.

439
440 All system suitability testing performed during the course of this study and all system suitability
441 test results thereof shall be reported to, and reviewed by, the PODP study coordinator before any
442 analytical data is accepted by the PODP Working Group. Failure to meet acceptance criteria will
443 be the basis for rejecting analytical data provided by the participating laboratory and frequent
444 failures by a participating laboratory can be the basis for the disqualification of that laboratory.

445
446 **Table 2. Composition of the System Suitability Test Mixtures.**

447
448 *Compounds for HPLC Analysis:*

449
450 Custom-made test mixture to be prepared by the participating laboratories from standard grade
451 reference materials:

452

Compound	Abbreviation	LC Test Mixture Concentration (µg/ml, ppm)
Caprolactam	CAP	1
Butylatedhydroxytoluene	BHT	5
Diphenylamine	DPA	5
Mono-(2-ethylhexyl) phthalate	MEHP	1
Stearic acid	SA	5
Di-(2-ethylhexyl phthalate)	DEHP	1
Bisphenol A	BPA	1

453

454 The test mix should be prepared by appropriate dilution of more concentrated stock solutions,
455 prepared using solvents appropriate for the individual reagents. The final composition of the test
456 mixture should be similar to, or compatible with, the mobile phase used in the LC analysis.
457

458 **Table 2. Composition of the System Suitability Test Mixtures (continued).**

459 *Compounds for GC Analysis, Grob Mixture:*

461 Commercial Sources:

462 e. g.: "Grob-Test-Mix", Cat# 11373, Restek

463 Reference:

464 K. Grob, Jr., G. Grob and K. Grob, "Testing Capillary Gas Chromatographic Columns", Journal
465 of Chromatography, 219, p. 13-20, (1981)
466
467
468

Combined solution of the following substances in methylene chloride:	Concentration, µg/ml (ppm)
	GC Test Mixture: (Grob Mixture diluted 1/20 in methylene chloride)
L(+)-2,3-butanediol	27
n-decane	14
2,6-dimethylaniline	16
2,6-dimethylphenol	16
methyl decanoate (C10:0)	21
methyl docecanoate (C12:0)	21
methyl undecanoate (C11:0)	21
nonanal	20
1-octanal	18
n-undecane (C11)	14

469 *Compounds for Headspace GC Analysis:*

470 Custom-made test mixture to be prepared by the participating laboratories from standard grade
471 reference materials:
472
473
474
475
476
477
478
479
480
481
482
483
484

Combined solution of the following substances in polyethylene glycol 200 ¹ (PEG 200):	HSGC Test Mixture I	
	µg/ml	µg/vial
Methanol	200	2
Acetic Acid	200	2
Cyclohexanone	100	1
Toluene	100	1
Trimethylsilanol ²	200	2
2-Ethyl hexanol	200	2

485 ¹ Preparation of SST-Sample:
486 - add 10 µl of the HS-Test-Mixture-I to a 20 ml crimp-cap vial
487 - add 10 µl of internal standard solution (2 mg of 1,4-Dioxane/ml PEG 200)
488 ²The material used is actually the sodium salt (sodium trimethylsilanolate).
489

490 The test mixture for headspace analysis can be prepared to contain the internal standard (1,4-
491 Dioxane) at the discretion of the testing laboratory.

492
493 *Composition of the ICP Test Mixture:*

494
495 System suitability testing for the ICP trace element analysis shall include the preparation and
496 testing of a system suitability test mixture that contains all the targeted elements listed previously
497 at a concentration of 0.25 mg/L.

498
499 The system suitability mixtures are minimally analyzed twice in the analytical runs, at the
500 beginning and at the end, thus establishing that adequate system performance is achieved and
501 maintained.

502
503 The evaluation of the system suitability results is as follows:

504
505 *LC Analysis:* The chromatograms for the system suitability test mixture are examined for the
506 presence of peaks corresponding to each analyte in the mix. While all analytes may not produce
507 responses in all detection methods, all analytes should produce peaks in at least one detection
508 method. All peaks should have a response with a signal to noise ratio (S/N) of 10 or greater.
509 The closest elution peak pair shall exhibit a resolution of greater than 1.5. All peaks should be
510 well-shaped, with a tailing factor less than 2.0. There should be no significant differences in the
511 chromatograms obtained at the beginning and the end of the chromatographic run. See Figure 2
512 for a sample chromatogram of the suitability test mixture.

513
514 *GC Analysis:* The chromatograms for the system suitability test mixture are examined for the
515 presence of peaks corresponding to each analyte in the mix. While all analytes may not produce
516 responses in all sample work-up methods (derivatized and non-derivatized), all analytes should
517 produce peaks in at least one work-up method. All peaks should have a response with a signal to
518 noise ratio (S/N) of 10 or greater. The closest elution peak pair shall exhibit a resolution of
519 greater than 1.5. All peaks should be well-shaped, with a tailing factor less than 2.0. There

520 should be no significant differences in the chromatograms obtained at the beginning and the end
521 of the chromatographic run. See Figure 3 for a sample chromatogram of the suitability test
522 mixture.


523
524 *HS-GC Analysis:* The chromatograms for the system suitability test mixture are examined for
525 the presence of peaks corresponding to each analyte in the mix. All analytes should produce
526 peaks that have a response with a signal to noise ratio (S/N) of 10 or greater. The closest elution
527 peak pair shall exhibit a resolution of greater than 1.5. All peaks should be well-shaped, with a
528 tailing factor less than 2.0. There should be no significant differences in the chromatograms
529 obtained at the beginning and the end of the chromatographic run. See Figure 4 for a sample
530 chromatogram of the suitability test mixture.

531
532
533
534
535
536
537
538
539
540
541
542
543
544
545
546
547
548
549
550
551
552
553
554
555
556
557
558
559
560
561
562
563
564

565
566
567
568
569
570

Figure 2. LC/UV/MS Chromatograms of the Suitability Mixture.


CAP = caprolactam; BPA = Bisphenol A; MEHP = mono-(ethylhexyl) phthalate; SA = stearic acid; DA = dehydroabiatic acid; DEHP = di-(2-ethylhexyl) phthalate. Peaks for BHT and DPA were not obtained in this run.



571

572 **Figure 3. GC/FID Chromatograms of the Grob Mixture.**


573
574 A. Underivatized.



575

Peak ID	Compound	Peak ID	Compound
1	2,3-Butanediol	7	2-ethyl hexanoic acid
2	Decane	8	2,6-Dimethyl aniline
3	1-Octanol	9	Methyl decanoate
4	Undecane	10	Dicyclohexylamine
5	1- Nonanal	11	Methyl undecanoate
6	2,6-Dimethyl phenol	12	Methyl dodecanoate

576 B. Derivatized.




577


Peak ID	Compound	Peak ID	Compound
1	Decane	7	2,6-Dimethyl phenol [TMS]
2	2,3-Butanediol [2TMS]	8	2,6-Dimethyl aniline [TMS]
3	Undecane	9	Methyl decanoate
4	1-Nonanal	10	Dicyclohexylamine
5	2-Ethyl hexanoic acid [TMS]	11	Methyl undecanoate
6	1-Octanol [TMS]	12	Methyl dodecanoate

Figure 4. GC/MS Chromatograms of the Headspace Suitability Mix.

578
 579
 580 Entire Chromatogram



581
 582
 583 Expanded view of 37 – 40 minutes.



584
 585

Retention Time (min)	Compound
16.1	Methanol
31.3	Trimethylsilanol
32.0	Toluene
33.2	1,4-dioxane
38.0	Cyclohexanone
39.1	Acetic acid
39.2	2-ethyl-1-hexanol

586

587 *ICP Analysis:* It shall be demonstrated that all elements can be detected at the prepared
588 concentration.
589

590 The performance expectations enumerated previously are general guidelines. All system
591 suitability data shall be reviewed by the Protocol's Study Coordinator and it is the responsibility
592 of the Coordinator to evaluate the system suitability data and establish its acceptability.
593

594 **C. Gas Chromatography (GC)**

595 *1. General*

596
597
598 Relatively volatile and semi-volatile compounds will be analyzed by Gas Chromatography (GC)
599 using a predominantly non-polar capillary column with wide (40 °C to 300 °C) temperature
600 programming. As noted previously, appropriate detection strategies will be employed (e.g. FID,
601 MS). Each GC analysis will produce an extractables "profile" in the form of a Total Response
602 Chromatogram (e.g. TIC for MS detection). As a first pass, identifications of individual
603 extractables will be accomplished with manual interpretation of the Electron Ionization (EI)
604 spectra assisted by computerized mass spectral library searching. Beyond this, more difficult
605 identifications may require the collection of additional data (such as Chemical Ionization GC/MS
606 for molecular weight confirmation and High Resolution Mass Spectrometry for elemental
607 composition), the purchase of reference compounds, *etc.* The PODP study coordinator shall be
608 consulted before a participating laboratory pursues the more difficult identifications.
609

610 *2. Sample Preparation*

611
612 The resulting extracts will usually contain low-level amounts of extractables. Sample
613 concentration and/or solvent switching may be necessary to provide compatible samples for the
614 analytical instrumentation. While it is possible to manipulate extracts to provide very large
615 concentration ratios, this has the undesirable effect of concentrating normal solvent impurities.
616 Therefore, extracts will be concentrated no more than 100X, which is reasonable given normal
617 ACS reagent purities of 99+%. The process for preparing (working-up) the aqueous extracts for
618 GC analyses is shown in Table 3. Similar evaporative sample concentration strategies may be
619 utilized with the organic extracts.
620
621
622
623
624
625
626
627
628
629
630
631

632

Table 3. Sample Work-up for Aqueous Extracts, GC Analysis	
Sample Preparation, Liquid-liquid Extraction; pH 2.5 and pH 9.5 Solutions.	<ol style="list-style-type: none"> 1 A 50-mL portion of each of the solutions is transferred to a 125 mL separatory funnel. 2 A 1.0-mL aliquot of the surrogate internal standard solution is added to each sample. 3 25 mL of Dichloromethane (DCM) is added to each funnel. 4 Each funnel is shaken for 1 minute. 5 The layers are allowed to separate and the lower (DCM) layer is collected. 6 Steps 3 through 5 are repeated. The collected DCM layers are combined. 7 The pH of each pH 2.5 sample is adjusted to ≈ 10 with 5 N NaOH. The pH of the pH 9.5 sample is adjusted to ≈ 2 with 5 N HCl. 8 Steps 3 through 5 are repeated twice for the pH adjusted samples. The collected DCM layers from all extractions are combined. 9 The DCM extracts are dried by adding anhydrous sodium sulfate to each collection flask. 10 Each DCM extract is transferred from the collection flask to a different Turbovap concentration tube with DCM rinses, and concentrated to less than 0.5 mL. A 0.5 mL aliquot of the injection internal standard is then added to the Turbovap tube. The final volume is adjusted to approximately 1 mL with DCM. 11 0.5 mL of each concentrated extract is transferred from the Turbovap tube to an autosampler vial. 12 The remaining 0.5 mL aliquot of each of dichloromethane extract described above is transferred to separate amber autosampler vials for TMS derivatization (see below)
Sample Preparation, Liquid-liquid Extraction; IPA/Water Solutions	The same basic process as noted above will be followed for the IPA/water samples. In the first extraction step, these samples will be pH adjusted to \approx pH 2 and extracted twice. In the second extraction step, the samples will be adjusted to \approx pH 10 and extracted twice. The resultant DCM extracts will be combined, dried and concentrated per steps 9 through 11 above.
TMS Derivatization of Residues	<ol style="list-style-type: none"> 1 Approximately 100 μL dimethyl formamide is added to each amber autosampler vial prepared under step 12 above. 2 The contents of each vial are evaporated nearly to dryness using nitrogen. 3 To each of the sample extracts, and the standard solutions is added 100 μL of BSTFA w/ 1% TMCS (Pierce) 4 Each vial is capped and allowed to stand for one hour at approximately 70°C. 5 DCM is added to each auto-sampler vial to make a final volume of approximately 0.5 mL, and is mixed.

633
 634
 635
 636
 637
 638
 639
 640
 641
 642
 643

The procedure contained in this Table is an example only and it is not required that participating laboratories adopt this procedure in either whole or in parts. However, any and all sample preparation procedures that will be used by a participating laboratory must be discussed with the PODP study coordinator prior to their utilization so that appropriate testing methodologies are utilized and harmonization between laboratories working on the same test articles can be achieved.

The procedure calls for the addition of a surrogate and injection internal standard, consistent with the system suitability assessment strategy enumerated previously. A surrogate internal standard is used to monitor the performance of the total procedure and is added to each extract in the initial

644 stage of its work-up. Requirements for such an internal standard are:


- 645
646 - sufficiently stable
647 - sufficiently soluble in all extraction solvents
648 - amenable to back-extraction from aqueous extracts by organic solvents
649 - semi-volatile
650 - amenable to all detection principles
651 - selectively detectable
652 - amenable to TMS-derivatization

653
654 The surrogate internal standard compound that meets these criteria has been identified as 4,4'-(m-
655 Phenylenediisopropylidene)diphenol (Bisphenol M):

656

CAS-no.: 13595-15-0
Molecular weight: 346.46
Molecular formula: $C_{22}H_{20}O_2$

Structure:



Source: e. g. Aldrich #450464

657
658 The Surrogate Standard Solution is prepared by dissolving 100 mg of Bisphenol M in 100 ml of
659 methanol, resulting in a concentration of 1000 µg/ml. This stock is further diluted 1 to 20 with
660 methanol to produce the surrogate internal standard solution containing 50 µg/mL Bisphenol M.


661
662 An injection internal standard is used to monitor the performance of the instrumental process
663 only and is added to each sample at the last stage of its work-up. Such an internal standard must
664 be:

- 665
666 - sufficiently stable
667 - sufficiently soluble in final extract
668 - semi-volatile
669 - amenable to all detection principles
670 - selectively detectable

671
672 The injection internal standard compound that meets these criteria has been identified as 4,4'-(m-
673 4,4'-Thiobis(3-methyl-6-t-butylphenol), Irganox 415:

CAS-no.: 96-69-5
Molecular weight: 358.538
Molecular formula: $C_{22}H_{30}O_2S$

Structure:



Source: e. g. Aldrich #366285

674
675 The Injection Standard Solution is prepared as follows: 100 mg of Irganox 415 are dissolved in
676 20 ml of methanol, concentration = 5000 µg/ml. This stock is further diluted 1 to 100 with
677 methanol to produce the surrogate internal standard solution containing 50 µg/mL Irganox 415.
678

679 The surrogate and injection internal standards are added to all samples to ensure that they are
680 properly worked-up and injected. Two internal standards are used to isolate the analytical
681 processes of sample work-up and instrumental analysis. The minimum performance expectation
682 for the internal standards is that they be present in the sample chromatograms with a response
683 whose signal to noise ratio is 10 or greater.
684

685 3. Operating Conditions

686
687 The following GC conditions (Table 4) serve as an illustration of a methodology which is
688 suitable for testing the prepared samples. The procedure contained in this Table is an example
689 only and it is not required that participating laboratories adopt this procedure in either whole or
690 in parts. However, any and all sample analysis procedures that will be used by a participating
691 laboratory must be discussed with the PODP study coordinator prior to their utilization so that
692 appropriate testing methodologies are utilized and harmonization between laboratories working
693 on the same test articles can be achieved.
694

695 Data cannot be collected while the injection solvent is in the ion source.
696

Table 4. Example Operating Parameters, GC Analysis of the Extracts.

Operating Parameter	Operating Value
Column	J&W DB-5HT, 30m x 0.25mm, 0.1 µm film thickness
Oven Program	Start at 50°C, hold for 5 min; ramp at 10°C/min to 300°C, hold for 5 min
Carrier Gas	He at 1 mL/min
Injection	Splitless; 2 µL.
Injector Temperature	310°C
FID Detector Temperature:	250°C
MS Transfer Line Temperature	310°C
MS Detection Details	70 eV (+), mass range of 33 – 650 amu (3.0 min or 6.0 min solvent delay used for un-derivatized or derivatized samples)

697

698 4. *General Comments.*

699
700 Note that the presented GC sample preparation and/or instrumental conditions are target
701 conditions for all participating laboratories and investigators. The actual conditions employed by
702 any participating laboratory should be reviewed by the PODP study coordinator prior to their
703 utilization so that appropriate testing methodologies are utilized and harmonization between
704 laboratories working on the same test articles can be achieved. In any event, the analyses
705 performed by the participating laboratory must meet system suitability criteria, as established in
706 Section VIII.A.

707
708 Any additional identification work beyond the first pass analysis will be performed only after
709 consultation with the PODP study coordinator.

710
711 Chromatograms of the extracts should be compared to chromatograms of the extraction blanks so
712 that peaks due to extractables can be delineated from peaks that reflect analytical artifacts.

713
714 **D. High Performance Liquid Chromatography (HPLC)**

715
716 1. *General*

717
718 Extracts and extraction blanks will be analyzed by High Performance Liquid Chromatography
719 with appropriate detection strategies, including DAD and MS as noted previously. The method
720 will use reversed-phase chromatography with a wide (gradient) range of solvent strengths. Each
721 LC analysis will produce several extractables “profiles” in the form of a Total Ion
722 Chromatogram (TIC), Extracted Ion Chromatograms (EIC) and UV chromatograms (total
723 response and/or specific UV wavelengths). As a first pass, identifications of individual
724 extractables will be accomplished with manual interpretation of the Atmospheric Pressure
725 Ionization Electrospray (API-ES) information. The LC and GC chromatograms will be
726 correlated to facilitate compound identification.

727
728 2. *Sample Preparation*

729
730 Unlike the GC analysis, the extracts and extraction blanks will typically not require extensive
731 sample preparation prior to HPLC analysis as the extraction matrices are generally compatible
732 with common HPLC mobile phases, thereby eliminating the need for solvent switching, and the
733 detection methods are sufficiently sensitive that sample concentration is not required. However,
734 some “solvent switching” may be necessary to produce samples that are HPLC-compatible.

735
736 3. *Operating Conditions*

737
738 The LC conditions in Table 5 serve as an illustration of a methodology which is suitable for
739 testing the prepared samples. The procedure contained in this Table is an example only and it is
740 not required that participating laboratories adopt this procedure in either whole or in parts.
741 However, any and all sample analysis procedures that will be used by a participating laboratory
742 must be discussed with the PODP study coordinator prior to their utilization so that appropriate

743 testing methodologies are utilized and harmonization between laboratories working on the same
744 test articles can be achieved.

745
746

Operating Parameter	Operating value	
Column	Agilent Zorbax Eclipse Plus C ₁₈ , 100 x 3.0 mm, 3.5µm particles	
Column Temperature	40°C	
Mobile Stage Components	A = 10 mM ammonium acetate, B = acetonitrile	
Mobile Stage Gradient	Time	% B
	0.0	5.0
	8.0	95.0
	11.0	95.0
	14.0	5.0
	17.0	5.0
Mobile Stage Flow Rate	0.8 mL/min	
Sample Size	60 µL	
Detection, UV	205 –300 nm	
Detection, MS	API-ES, positive ion and negative ion (mass range 80 – 1200)	
Sample Preparation	None, direct injection	

747

748 4. *General Comments*

749

750 Any additional identification work beyond the first pass analysis will be performed only after consultation
751 with the PODP study coordinator.

752

753 Chromatograms of the extracts should be compared to chromatograms of the extraction blanks so
754 that peaks due to extractables can be delineated from peaks that reflect analytical artifacts.

755

756 **E. Headspace GC/MS (HS-GC)**

757

758 1. *General*

759

760 Headspace analysis of extracts allows for an assessment of the volatile organic extractables.
761 Volatiles present in the extract are thermally evolved into the headspace. The evolved volatile
762 entities are “captured” in the headspace gas, which is transferred, in whole or in part, to an
763 appropriate analytical technique. Since the headspace sample is a gas, gas chromatography is the
764 analytical method of choice. Mass spectrometry is the detection method of choice because it
765 facilitates the identification of evolved entities.

766

767 The headspace methodology is intended to uncover volatile entities that are present in the
768 extract; it is not intended to produce “volatiles” by causing extractables present in the extracts to
769 thermally decompose. Thus the headspace “extraction” is accomplished at relatively low
770 temperatures (e.g. 120°C or less).

771

772

773 2. *Sample Preparation*

774
775 Place approximately 4 mL of sample (extract or extraction blank) into a 20 mL headspace
776 autosampler vial containing approximately 10 grams anhydrous sodium sulfate. Seal the vial by
777 crimping a cap onto it.

778
779 For semiquantitative evaluation and also to check for proper performance of the measurement, an
780 10 µL aliquot of a solution of 1,4-Dioxane in polyethylene glycol 200 (concentration 2 mg/mL)
781 is added to each vial. This solution is prepared as follows: 20 mg of 1,4-Dioxane are dissolved in
782 10 mL of polyethylene glycol 200, resulting in a concentration of 2 mg/mL.

783
784 Note: A positive displacement pipetting system (e. g. Gilson Microman[®]) should be used for
785 dosing this solution due to its high viscosity.

786
787 3. *Operating Conditions*

788
789 The operating conditions for the Headspace GC/MS are contained in Table 6.
790

Table 6. Operating Parameters, Headspace GC/MS Analysis for Volatiles.	
Operating Parameter	Operating Value
A. Headspace Autosampler	
Oven Temperature	80°C
Needle Temperature	120°C
Transfer Line Temperature	155°C
Carrier gas	He at 5 psi
Equilibrium Time	120 min
B. GC/MS Analyzer	
Column	J&W DB-WAXETR, 60 m x 0.32 mm I.D., 1 µm film
Oven Program	Start at 35°C, hold for 7 minutes. Ramp at 1°C/min to 40°C, hold for 15 minutes. Ramp at 10°C/min to 100°C. Ramp at 25°C/min to 240 °C, hold for 5 min.
MS Ionization Mode	EI+, 70 eV
MS Transfer Line Temperature	240°C
MS Detection Mass Range	25 – 200 amu
Solvent Delay	0 min

791
792 4. *General Comments*

793
794 The analyses performed by the participating laboratory must meet system suitability criteria, as
795 established in Section VIII.A.

796
797 The Headspace GC/MS analysis will produce an extractables “profile” in the form of a Total
798 Response Chromatogram (e.g. TIC for MS detection). As a first pass, identifications of
799 individual extractables will be accomplished with manual interpretation of the Electron

800 Ionization (EI) spectra assisted by computerized mass spectral library searching. More difficult
801 identifications may require the collection of additional data (such as Chemical Ionization GC/MS
802 for molecular weight confirmation and High Resolution Mass Spectrometry for elemental
803 composition), should be discussed with the PODP study coordinator before a participating
804 laboratory pursues these more difficult identifications.

805
806 Chromatograms of the extracts should be compared to chromatograms of the extraction blanks
807 (Headspace vials containing no test material) so that peaks due to extractables can be delineated
808 from peaks that reflect analytical artifacts.

809
810 The concentration of any extractables can be estimated via the use of the internal standard.

811
812 **F. Inductively Coupled Plasma Atomic Spectroscopy (ICPAS)**

813
814 *1. General*

815
816 Single elements (e.g. metals) in the extracts will be analyzed by Inductively Coupled Plasma
817 Atomic Spectroscopy using appropriate methods and techniques for the determination of
818 common analytes. Detection strategies such as optical emission and mass spectrometry shall be
819 employed. ICP analyses should be performed consistent with USP practices.⁴

820
821 *2. Sample Preparation*

822
823 The resulting extracts will usually contain low-level amounts of extractables as well as the ionic
824 constituents of the extracting media. The pH 2.5 and pH 9.5 extracts will contain large quantities
825 of sodium and the pH 9.5 extracts will have large quantities of phosphorous. Thus these analytes
826 cannot be determinable in these extracts.

827
828 The material extracts will need to be processed to some extent prior to testing. Aqueous samples
829 will be acidified directly via addition of nitric acid. Place 10 mL of aqueous extract in a trace
830 metal-free plastic vessel. Add 0.5 mL of concentrated nitric acid and mix well.

831
832 *3. Operating Conditions*

833
834 The ICP spectrometers shall be operated consistent with good laboratory practices and standard
835 procedures in place in the participating testing facilities. The following is a list of elements that
836 must be included in the ICP analysis: Al, As, Be, B, Cd, Cr, Co, Cu, Fe, Li, Mg, Mn, Mo, Ni, Pd,
837 Pt, S, Sb, Se, Si, Sn, Sr, Ti, V, W, Zn, and Zr. Additional elements may be reported depending
838 on the capabilities of the lab. The analysis conditions should be such that these elements can be
839 measured at the appropriately low levels, typically 0.01 µg/mL or less in the material extracts.

840
841
842
843

⁴ USP 30, <730> Plasma Spectroscopy.

844 4. *General Comments*

845
846 The analyses performed by the participating laboratory must meet system suitability criteria, as
847 established in Section VIII.B.

848
849 Results for the extracts should be compared the results for the extraction blanks so that
850 extractables can be delineated from analytical artifacts and solvent impurities/components.

851

852 **IX. DATA EVALUATION AND REPORTING**

853

854 **A. Qualitative Analysis**

855

856 • A list of all identified entities (compounds, elements) that were not detected in the
857 corresponding blank. This list should include the recognized compound name,
858 CAS Registry number, chemical formula, and chemical structure.

859 • A list of all unidentified chromatographic peaks that were not detected in the
860 corresponding blank at signal to noise ratios greater than 10. The participating
861 laboratory should determine and report the analyte concentration that corresponds
862 to this signal to noise ratio (typically defined as the limit of quantitation, LOQ).

863 • Copies of chromatograms, spectra, etc.

864 • Complete methodological information for both the extraction and analysis
865 processes.

866 • The required system suitability results, which should include an assessment of
867 detectability.

868 • The identification status for all compounds shall be established and reported as
869 follows:

870

871 • A *Confirmed* identification means that collaborating information has been
872 obtained including mass spectrometric fragmentation pattern, confirmation of
873 molecular weight (or elemental composition), match in retention time and
874 spectrum with authentic standard.

875 • A *Confident* identification means that sufficient data to preclude all but the most
876 closely related structures have been obtained

877 • A *Tentative* identification means that data have been obtained that are consistent
878 with a class of molecule only.

879

880 • A report format will be distributed to the participating laboratories.

881

882 **B. Semi-Quantitative Analysis**

883

884 While it is not the primary intent of this Stage 1 Protocol to produce quantitative data,
885 some of the test methods employed may be amenable to concentration estimation (e.g.
886 ICP, GC with internal standards). In the case that a participating laboratory reports
887 concentration estimates, the means by which such estimates were obtained must be
888 indicated. Additionally, all such estimates shall be reported with a convention (e.g.

889 significant figures) which effectively reflects the uncertainty in the determination. As
890 was noted previously, the threshold for reporting semi-quantitative results is 0.1 ppm
891 ($\mu\text{g/ml}$).
892

893 X. GLOSSARY OF ABBREVIATIONS

	GC/FID	Gas Chromatography with Flame Ionization Detector
	GC/MS	Gas Chromatography with Mass Spectrometric Detection
	HPLC/DAD	High Pressure Liquid Chromatography-Diode Array Detection
	LC/MS	Liquid Chromatography Mass Spectrometric Detection
	ICP/AES	Inductively Coupled Plasma Atomic Emission Spectroscopy
	ICP-MS	Inductively Coupled Plasma Mass Spectrometry
	HS-GC	Gas Chromatography with Headspace gas Sampling
	TIC	Total Ion Chromatogram
	API-ES	Atmospheric Pressure Ionization - Electrospray
	HS	Headspace
894	PQRI	Product Quality Research Institute
895	OINDP	Orally Inhaled and Nasal Drug Products
896	PODP	Parenteral and Ophthalmic Drug Products
897	BFS	Blow-fill-seal
898		

899 XI. REFERENCES

900
901 PQRI Research Project Proposal: *Development of Scientifically Justifiable Thresholds and Best*
902 *Demonstrated Characterization Practices for Leachables and Extractables in Parenterals and*
903 *Ophthalmic Drug Products (PODP)*. (March 2008).

904
905 FDA Guidance for Industry: *Container Closure Systems for Packaging Human Drugs and*
906 *Biologics-Chemistry Manufacturing, and Controls Documentation*, (May 1999).

907
908 European Medicines Agency (EMA): *Guideline on Plastic Immediate Packaging Materials*,
909 (May, 2005).

910
911 **Parenteral and Ophthalmic Drug Products Leachables and Extractables Working Group.**
912 *Experimental Protocol for Qualitative Controlled Extraction Studies on Material Test Articles*
913 *Representative of Prefilled Syringe (PFS) and Small Volume Parenteral (SVP) Container*
914 *Closure Systems*. Study Protocol – Stage 1. Issued and Effective, December, 2009; Amendment
915 #1, September, 2011.
916